

BUNION PROTECTOR with EPITHELIUM 26®

INDICATION

Using the EPITACT® bunion protector relieves the pain associated with the pressures and frictions on the bunion.

CONTRA-INDICATIONS

Do not use on injured skin. In case of arthritis, diabetes, lower limb neuropathy or ulcers, seek advice from your doctor or podiatrist before use and then regularly monitor your foot.

PRECAUTIONS FOR USE

- The hallux valgus protection pad should be worn next to the skin, on the right or left foot.
- Slip your foot into the protection pad, putting your great toe in the tubular part. The Epithelium 26® gel module should be placed on the hallux valgus.
- To help the protection pad to remain in position, wear tights, socks or stockings over it.

FEATURES

For bunion protection, size matters: carrying a protection that is too bulky in the shoe can aggravate the problem further or create undesirable pain on other areas of the foot. The EPITACT® protector is therefore thin and lightweight, but still very effective. The patented EPITHELIUM 26® silicone

gel helps avoid frictions on the bunion and was deliberately limited to 1mm of thickness.

DESIGN

Thin and lightweight, the protector can be worn in any shoe.

ADVICE

If your great toe is still mobile laterally, we recommend wearing the EPITACT® flexible bunion corrector. It limits the evolution of the deformation, relieves joint pain and avoids rubbing/friction on the bunion.

CLEANING & FIFESPAN

Wash in a machine at 30°C / 86°F or by hand with soap and water (without rubbing) and leave to dry naturally. Do not iron. Reusable for several months.

COMPOSITION: 61% polyamide/elastane fabric, 39% silicone.

The materials used have passed all the regulatory tests required by ISO standard 10993-1 relating to products that come into contact with the skin. If an incident occurs, stop using the product immediately and inform MILLET Innovation. Millet Innovation cannot be held responsible for failure to use the product in accordance with the manufacturer's recommendations. French is the reference language. CE initial marking: 2001.



Machine washable
at 30°C / 86°F.



Do not iron
or steam.



Do not dry clean.



Do not bleach.



Do not
tumble-dry.



See packaging.

Product of MILLET Innovation research.



MILLET Innovation ZA Champgrand
26270 Loriol/Drôme - FRANCE

MILLET
INNOVATION

www.epitact.co.uk