

## HAMMER TOE CUSHIONS with EPITHELIUM26®

### INDICATIONS


Hammer, claw toes, corns at the end of the toes, high-arched feet.

### CONTRA-INDICATIONS

Do not use on injured skin. In case of arteritis, diabetes, lower limb neuropathy or ulcers, seek advice from your doctor or podiatrist before use and then regularly monitor your foot. The strap may cause friction between the two toes. This solution should therefore be avoided by diabetics suffering from neuropathy.

### TIPS FOR USE

- Hold the hammer toe cushion by the rounded top part. The elastic strap slips over the third toe. The largest part of the cushion needs to be under the second toe.
- Place it in the free space at the base of the toes.
- Then slip on socks, stocking or tights to secure the cushion in proper position.



### FEATURES

The EPITACT® Hammer Toes Cushions straighten and lift the toes softly to relieve pressures and frictions at the tips.

### DESIGN

The hammer toe cushions are filled with a material, covered by soft and breathable fabric, called EPITHELIUM™.

### CLEANING & LIFESPAN

Wash in a machine at 30°C / 86°F or by hand with soap and water (without rubbing) and leave to dry naturally. Do not iron. Reusable for several months.

### COMPOSITION: 90% silicone, 10% polyamide/elastane fabric.

The materials used have passed all the regulatory tests required by ISO standard 10993-1 relating to products that come into contact with the skin. If an incident occurs, stop using the product immediately and inform MILLET Innovation. Millet Innovation cannot be held responsible for failure to use the product in accordance with the manufacturer's recommendations. French is the reference language. CE initial marking: 2003.



Machine washable at 30°C / 86°F.



Do not bleach.



Do not iron or steam.



Do not tumble-dry.



Do not dry clean.



See packaging.

Product of MILLET Innovation research.



MILLET Innovation ZA Champgrand  
26270 Loriol/Drôme - FRANCE

MILLET  
INNOVATION

www.epitact.co.uk

CE

NT0634EN-J  
09/2018